

When One Defibrillating Shock Isn't Enough

New Clinical Evidence for Increasing the Dosage


Escalating dosage: a strategy to maximize defibrillation effectiveness

New clinical evidence from two out-of-hospital ventricular fibrillation (VF) studies^{1,2} confirms the effectiveness of current AHA recommendations for first shock dosage. But for patients who need additional shocks, these studies show that repeating the same first shock dosage is inferior to a strategy of increasing to a higher dosage.



An escalating protocol is consistent with the AHA Guidelines 2005 recommendation for subsequent dosage;

“...it is reasonable to use selected energies of 150J to 200J with a biphasic truncated exponential waveform...for the initial shock. For second and subsequent shocks use the same or higher energy (Class IIa)³.” (emphasis added)

Furthermore, this new clinical evidence from these two studies shows a benefit from escalating dosage.



- High rate of 1st shock success with LIFEPAK® Defibrillator/Monitor at (92%)¹
- Diminishing return from repeating the dosage after a first shock failed (92% first shock vs. 61% second shock success, P=0.001)
- 121/467 (26%) known survivors at 30 days: 51% of them received 360J



- Triple blinded, randomized controlled clinical trial in 221 patients comparing two biphasic dosage protocols: fixed lower energy 150J vs. escalating higher energy 200J-300J-360J.²
- Among patients requiring more than one shock, the escalating higher energy regimen provided significantly higher rates of conversion to an organized rhythm (36.6% vs. 24.7%, P=0.035) and VF termination (82.5% vs. 71.2%, P=0.027) compared to a fixed low energy regimen.

When One Defibrillating Shock Isn't Enough

New Clinical Evidence for Increasing the Dosage

Broad Dosage Capability: An Important Consideration for Patient Care

Clinical studies have shown a majority of cardiac arrest victims with an initial rhythm of VF will experience repeated episodes of VF over the course of a resuscitation attempt.⁴⁻⁸ For any given VF episode the initial shock may fail; in fact defibrillation becomes more difficult with each recurrent episode.¹ For these patients in both hospital and out-of-hospital settings, increasing the dosage of subsequent shocks above the dose used for the first shock has proven to be a better strategy for terminating VF than simply repeating a failed dosage.^{1,2,9} LIFEPAK defibrillator/monitors provide broad dosage capability, up to 360J. For more information, please contact your Physio-Control representative at 800.442.1142 or visit www.physio-control.com.

REFERENCES

- 1 Chapman FW, Walker RG, Koster RW. Use of 360 Joule Biphasic Shocks for Initial and Recurrent Ventricular Fibrillation in Prehospital Cardiac Arrest. *Resuscitation* 2006;69(1):49-50.
- 2 Stiell IG, et al. The BIPHASIC Trial: A randomized comparison of fixed lower versus escalating higher energy levels for defibrillation in out-of-hospital cardiac arrest. *Circulation* 2007 115: 1511-1517.
- 3 Guidelines 2005 for Cardiopulmonary Resuscitation and Emergency Cardiovascular Care. *Circulation* 2005; 112(24): IV-37.
- 4 Stothert JC, et al. Rectilinear biphasic waveform defibrillation of out-of-hospital cardiac arrest. *Prehospital Emergency Care* 2004;8(4):388-92.
- 5 Gliner BE, et al. Treatment of out-of-hospital cardiac arrest with a low-energy impedance-compensating biphasic waveform automatic external defibrillator. *Biomedical Instrumentation and Technology* 1998;32:631-44.
- 6 van Alem AP, Post J, Koster RW. VF recurrence: characteristics and patient outcome in out-of-hospital cardiac arrest. *Resuscitation* 2003;59:181-88.
- 7 White RD, Russell JK. Refibrillation, resuscitation and survival in out-of-hospital sudden cardiac arrest victims treated with biphasic automated external defibrillators. *Resuscitation* 2002;55:17-23.
- 8 Weaver WD, et al. Ventricular defibrillation – a comparative trial using 175-J and 320-J shocks. *N Engl J Med* 1982;307:1101-6.
- 9 Walsh SJ, et al. Efficacy of distinct energy delivery protocols comparing two biphasic defibrillators for cardiac arrest. *Am J Cardiol* 2004;94:378-80.

For further information please contact your local Physio-Control representative or visit www.physio-control.com



HEADQUARTERS / MANUFACTURING

Physio-Control, Inc.
11811 Willows Road NE
P. O. Box 97006
Redmond, WA 98073-9706 USA
Tel 425 867 4000
Fax 425 867 4121
www.physio-control.com

SALES OFFICES

Physio-Control UK
Medtronic Ltd
Suite One, Sherbourne House
Croxley Business Park
Watford, Herts
WD18 8WW
Tel 44 1923 212 213
Fax 44 1923 241 004

Canada
Medtronic of Canada
6733 Kitimat Road
Mississauga, ON
L5N 1W3
Tel 888 879 0977
Fax 416 430 6115