

LUCAS® 3 Chest Compression System

Myth Busted! The LUCAS device does fit large patients!

Weight:

82 kg / 12.9 stone

Height:

183 cm



Weight:

150 kg / 23.6 stone

Height:

196 cm



Weight:

145 kg / 22.9 stone

Height:

178 cm



Proper Placement is on Upper Chest

Do not assume the LUCAS Chest Compression System can't fit larger patients. Many large patients carry their weight through their mid-section. When LUCAS is properly applied, it rests over the upper chest or sternum, with enough room for the compressor to give compressions according to guidelines. As long as the unit snaps into the backboard and the suction cup does not compress the patient's chest in the start position, it will operate as intended.



Patients Eligible for Treatment:

- Sternum height of 17 – 30.3 cm
- Maximum chest width: 45 cm

Using the LUCAS device is not restricted by patient weight.

LUCAS® 3 Chest Compression System

Myth Busted! The LUCAS device does fit large patients!



For further information please contact your local Physio-Control representative or visit our website at www.physio-control.com



Physio-Control Headquarters
11811 Willows Road NE
Redmond, WA 98052 USA
Tel 425 867 4000
Fax 425 867 4121
www.physio-control.com

**Physio-Control Operations
Netherlands B.V.**
HQ Europe, Middle East & Africa
UP Office Building
Piet Heinkade 55,
1019 GM Amsterdam NL
Tel +31 (0)20 7070560
Fax +31 (0)20 3391194
www.physio-control.com

**Physio-Control
UK Sales Ltd**
11 Old Jewry, 7th Floor
London EC2R 8DU
United Kingdom

**Physio-Control
Australia Pty Ltd**
Suite 4.01
15 Orion Road
Lane Cove
NSW 2066
Australia
Toll Free Tel 1800 987 982
Toll Free Fax 1800 890 892

 **Jolife AB**, Scheelevägen 17, Ideon Science Park, SE-223 70 LUND, Sweden