



LIFEPAK[®] 15

MONITOR / DEFIBRILLATOR

Genuine Accessories from Physio-Control

LIFEPAK[®] 15 MONITOR/DEFIBRILLATOR





Thank you for choosing Physio-Control as your partner in helping save lives and improve patient care.

You already know our reputation for quality and reliability in LIFEPAK defibrillators and AEDs. When it comes to LIFEPAK accessories, we know you won't accept anything less. Inside this catalog you'll find accessories designed to meet your needs and reliability standards. All the accessories and disposables in this catalog are designed to work with the LIFEPAK 15 monitor/defibrillator. If you have any questions about accessories or disposables, please contact your Physio-Control representative.

Temperature Monitoring



Temperature Adapter Cable

11140-000079 (10ft)
11140-000078 (5ft)



Esophageal-Rectal Temperature Sensor

11996-000360 (9FR, 20/box)



Skin Temperature Sensor

11996-000359 (20/box)



Foley Catheter Temperature Sensor

11996-000361 (1FR, 10/pk)
11996-000362 (16FR, 10/pk)
11996-000363 (18FR, 10/pk)

NIBP Monitoring Accessories



NIBP Tubing

21300-007299 (9ft)
21300-007298 (12ft)



Reusable Cuff

X-Large Adult
 45 - 44 cm
11160-000009

Pediatric
 13 - 20 cm
11160-000003



Single Patient Use Cuff

X-Large Adult
 45 - 44 cm
11160-000010

Pediatric
 13 - 20 cm
11160-000004



NIBP Tubing, Coiled

21300-007300 (2-9ft)

Large Adult
 32 - 42 cm
11160-000007

Infant
 8 - 14 cm
11160-000001

Large Adult
 32 - 42 cm
11160-000008

Infant
 8 - 14 cm
11160-000002

Adult
 26 - 35 cm
11160-000005

Adult
 26 - 35 cm
11160-000006

Power Options



REDI-CHARGE Base
AC power cord and adapter tray not included.
11141-000115

REDI-CHARGE LIFEPAK 15 Adapter Tray
11140-000052

AC Power Cord
11140-000015



Lithium-ion Battery 5.7 Ah
5.7 amp hour, 11.1 volt, rechargeable, with fuel gauge.
21330-001176



Mobile Battery Charger
Includes AC and DC power cords, mounting bracket and operating instructions.
11577-000011



AC Power Adapter
Includes Right Angle Cable (AC power cord not included)
11140-000072

DC Power Adapter
Includes DC cable and Right Angle Cable
11140-000074

AC Power Cord
11140-000015

Power Attachment Kit
11577-000019



Extension Cable for AC/DC Power Adapter
11140-000080



Replacement Right Angle Power Cable for AC/DC Power Adapter
11140-000081



Replacement DC Input Cable
11140-000084

ECG Monitoring Accessories



**12-Lead ECG Cable Trunk
Cable with 4-Wire Limb Leads**
11111-000018 (5ft)
11111-000020 (8ft)



**12-Lead ECG Cable 6-Wire
Precordial Attachment**
11111-000022




4-wire Cable Comb
21300-008054 (10/pk)

6-wire Cable Comb
21330-008055 (10/pk)



3-Wire ECG Cable
Right-angle connector
11110-000029



5-Wire ECG Cable
Right-angle connector, 4-wire limb plus one chest lead, labeled "V1" on the LIFEPAK 15 monitor reports.
11110-000066



Box of Strip Chart Recorder Paper
100mm x 22m
11240-000016 (2 rolls/box)



LIFE-PATCH[®] ECG Electrodes
Adult, pregelled.
11100-000001 (3/pk)



LIFE-PATCH ECG Electrodes
Adult, pregelled.
11100-000002 (4/pk)

Therapy Delivery Accessories

Hard Paddles and Electrode Gel



Standard Hard Paddles

11130-000061 (1 pair)



Pediatric Paddle, External

(Two required); slips onto standard adult hard paddle.

11133-000007 (1 paddle)



SIGNAGEL® Electrode Gel

For use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback and EMG.

21300-005847 (8.5 oz./tube)



Internal Paddles (Requires Internal Paddle Handles and Internal Paddles Adapter Cable)

1 in. size

11131-000010 (1 pair, 6.25 in. shaft)

1.5 in. size

11131-000011 (1 pair, 6 in. shaft)

11131-000021 (1 pair, 9 in. shaft)

11131-000024 (1 pair, 14 in. shaft)

2 in. size

11131-000012 (1 pair, 5.75 in. shaft)

11131-000022 (1 pair, 8.75 in. shaft)

2.5 in. size

11131-000013 (1 pair, 5.5 in. shaft)

11131-000019 (1 pair, 8.5 in. shaft)

3.5 in. size

11131-000014 (1 pair, 5 in. shaft)

11131-000023 (1 pair, 8 in. shaft)



Internal Paddle Handles with Discharge Control

(For use with the Internal Paddles Adapter Cable)

11131-000001 (1 pair)



Internal Paddles Adapter Cable

(For use with Internal Paddle Handles)

111998-000326

Therapy Delivery Accessories

EDGE System™ Electrodes for Pacing/Defibrillation/ECG with QUIK-COMBO® Connector

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector

24 in. leadwire length

11996-000091



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

24 in. leadwire length

11996-000090



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System

42 in. leadwire length

11996-000017



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators;
12 month minimum shelf life at time of shipment
24 in. leadwire length.

11996-000093



QUIK-COMBO Therapy Cable

With convenient TRUE-LOCK™ Cable Connector.
Length is approximately 8 ft. For use with
LIFEPAK 15 monitor/defibrillator.

11113-000004

Pulse Oximetry Monitoring Accessories

Masimo SET® RC Patient Cables



RC Patient Cable

For use with M-LNCS and Rainbow Patient Sensors

11171-000037 (4ft)

11171-000038 (12ft)

Masimo SET RC Patient Cable Compatible SpO₂ Sensors



M-LNCS Reusable Sensor

11171-000046 (Ad)

11171-000047 (Ped)



M-LNCS Adhesive Sensors (20/box)

11171-000039 (Ad)

11171-000040 (Ped)



M-LNCS Adhesive Sensors (20/box)

11171-000043 (Neo/Pt)

11171-000042 (Neo/Ad)

11171-00004A (Inf)

Masimo SET RC Patient Cable Compatible Rainbow® SpO₂, SpCO, SpMet Sensors



Rainbow Reusable Sensor

11171-000049 (Ad)

11171-000050 (Ped)



Rainbow Adhesive Sensor (10/box)

11996-000342 (Inf)

11996-000341 (Neo/Ad)



Rainbow Adhesive Sensor (10/box)

11996-000339 (Ad)

11996-000340 (Ped)

Pulse Oximetry Monitoring Accessories

Masimo SET LNC Patient Cables



Red LNC Patient Cable

For Use with LNCS Patient Sensors

11996-000323 (4ft)

11996-000324 (10ft)

11996-000325 (14ft)

Masimo SET LNC Patient Cable Compatible SpO₂ Sensors



LNCS® Reusable Sensor

11171-000017 (Ad)

11171-000018 (Ped)



LNCS Reusable Soft Sensor

11171-000052 (Ad)



LNCS Adhesive Sensor (20/box)

11171-000019 (Ad)

11171-000020 (Ped)



LNCS Adhesive Sensor (20/box)

11171-000029 (Neo/Pt)

11171-000028 (Neo/Ad)

11171-000031 (Inf)

Pulse Oximetry Monitoring Accessories

LIFEPAK 15 Direct Connect SpO₂ Only Patient Sensors



Adult Reusable Direct Connect Sensor

11996-000331 (3ft)
11996-000332 (12ft)



Pediatric Reusable Direct Connect Sensor

11996-000333 (3ft)
11996-000334 (12ft)



Adult Reusable Soft Direct Connect Sensor

11171-000053 (8ft)

LIFEPAK 15 Direct Connect Rainbow SpO₂, SpCO, SpMet Patient Sensors



Adult Rainbow Direct Connect Reusable Sensor

11996-000335 (3ft)
11171-000032 (8ft)
11996-000336 (12ft)



Pediatric Rainbow Direct Connect Reusable Sensor

11996-000337 (3ft)
11171-000033 (8ft)
11996-000338 (12ft)

Additional Masimo Accessories



Red MNC Cable

Connects LIFEPAK 15 to Nellcor patient sensor
11996-000365 (4ft)
11996-000366 (10ft)



Reusable Ambient Light Shield

11171-000054 (5/bag)



Disposable Ambient Light Shield

11171-000055 (10/bag)

Pulse Oximetry Monitoring Accessories

Nellcor Reusable Sensors and Cables



DURASENSOR Reusable Clip
11996-000060 (Adult - Ref# DS100A)

DURA-Y Multisite Reusable Sensor
11996-000106 (≥ 1 kg - Ref# D-YS)



Oxiband Reusable Sensor
Includes 50 disposable adhesive sensors
11996-000061 (Ad/Neo - Ref# OXI A/N)
11996-000062 (Ped/Inf - Ref# OXI P/I)

Disposable Adhesive Bandage Wrap
Not Available in Canada. (100/pk)
11996-000048 (Ad/Neo - Ref# ADH A/N)
11996-000049 (Ped/Inf - Ref# ADH P/I)

Disposable Adhesive Foam Wrap
Not Available in Canada. (100/pk)
11996-000110 (Ad/Neo - Ref# FOAM A/N)
11996-000108 (Ped/Inf - Ref# FOAM P/I)

Nellcor Disposable Sensors

Single patient use



Oxisensor II (24/box)
Adult
11996-000113 (18 in. - Ref# D25)
11996-000114 (36 in. - Ref# D25L)

Pediatric
11996-000116 (18 in. - Ref# D20)



Oxisensor II (24/box)
Infant
11996-000115 (18 in. - Ref# I20)

Neonatal/Adult
11996-000117 (18 in. - Ref# N25)

End-Tidal CO₂ (EtCO₂) Monitoring Accessories

Oridion® Filterlines for Intubated Patients

Single patient use



FilterLine® SET

Key Applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

11996-000081 (25/pk, 200 cm)

11996-000164 (25/pk, 400 cm)



FilterLine H SET

Key Applications: Critical care, Humidified Environments.

Adult/Pediatric

11996-000080 (25/pk, 200 cm)

Adult/Neonatal

11996-000001 (25/pk, 200 cm)

Oridion Non-Intubated Filterlines

Single patient use



Smart CapnoLine® Plus

Key Applications: Procedural Sedation, Upper GI Procedures, MAC, RMS, BD, Rapid Response Teams.

Adult with O₂

11996-000163 (25/pk, 200 cm)

11996-000167 (100/pk, 200 cm)

11996-000165 (25/pk, 400 cm)

Adult without O₂

11996-000162 (25/pk, 200 cm)

11996-000166 (100/pk, 200 cm)



Smart CapnoLine

Key Applications: Procedural Sedation, SMS, SD, Rapid Response Teams.

Pediatric with O₂

11996-000128 (25/pk, 200 cm)

Pediatric without O₂

11996-000120 (25/pk, 200 cm)

Cases & Mounting Options



Standard Carrying Case

Includes right pouch and left pouch, for use with LIFEPAK 15 monitor/defibrillator.
11577-000002

Top Pouch

Storage for sensors and electrodes; insert in place of standard paddles.
11220-000028

Shoulder Strap

For use with LIFEPAK 15 monitor/defibrillator.
11577-000001

Back Pouch

Ideal for additional accessory storage.
For use with LIFEPAK 15 monitor/defibrillator.
11260-000039



Bed Connector

For use in hospital only.
11996-000374

Communication Accessories



LIFEPAK Monitor to PC Cable

For connecting LIFEPAK 12 or LIFEPAK 15 monitor/defibrillator to PC.
11230-000020 (Serial)
11996-000369 (USB)



Titan II Wireless Gateway

For transmitting data from LIFEPAK 12/15 to the LIFENET System or CODE-STAT data review software. Requires existing wireless network.
21996-000073
21996-000075 (with voice recording)



3G Gateway

For transmitting data from LIFEPAK 12/15 to the LIFENET System or CODE-STAT data review software. Requires data plan.
21996-000086 (Verizon Gateway)
99428-000305 (Verizon Data Plan)
21996-000082 (AT&T Gateway)
99428-000304 (AT&T Data Plan)

Testers & Training Materials



12-Lead Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias, and bradycardias, as well as ST segment and T wave abnormalities. For use with LIFEPAK devices with a 12-lead ECG feature.

11996-000311

3-Lead Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe, interactive training. Select from 17 ECG rhythms including: fibrillation, tachycardias and bradycardias.

11996-000310



Defibrillator Checker

Tests integrity of energy delivery through standard hard paddles. Neon light indicates energy has been delivered.

11998-000060



Test Load

Use to perform therapy cable performance checks. Connects to the QUIK-COMBO therapy cable on the defibrillator.

21330-001365 (English only)

21330-001367 (English, French for Canada)

Training Tools



Clip-on Training Electrodes

For use with QUIK-COMBO patient simulators. Allows users to practice defibrillation electrode placement. Includes instructions for use and reusable training electrodes.

11250-000052 (5/pair)

Replacement Training Electrodes

Reusable training electrodes. Adhesive electrode portion only.

11101-000003 (5/pair)



Inservice Video: LIFEPAK 15 Monitor/Defibrillator

21330-001357 (DVD - NTSC version)

Operating Instructions: LIFEPAK 15 Monitor/Defibrillator

Also available as free online download at www.physio-control.com

26500-002408



UNIVERSITY

Physio-Control University

12-Leads Made Easy

44500-000001

Capno Made Easy

44500-000003

For further information, please contact Physio-Control at 800.442.1142 (U.S.), 800.895.5896 (Canada) or visit our website at www.physio-control.com.



Physio-Control Headquarters
Redmond, WA 98052
www.physio-control.com

Customer Support
P. O. Box 97006
Redmond, WA 98073
Toll free 800 442 1142
Fax 800 426 8049

Physio-Control Canada
Physio-Control Canada Sales, Ltd.
7111 Syntex Drive, 3rd Floor
Mississauga, ON
L5N 8C3
Canada
Toll free 800 895 5896
Fax 866 430 6115