



LIFEPAK[®] 20

DEFIBRILLATOR/MONITOR SERIES

Genuine Accessories from Physio-Control

LIFEPAK 20 MONITOR/DEFIBRILLATOR SERIES





Genuine Accessories from Physio-Control

Ensure the safety of your staff and patients. Only Physio-Control accessories and disposables have been thoroughly tested with LIFEPAK products. Stick with the name you trust and the company that stands behind it. To order, contact your Physio-Control representative.

Power Options



AC Power Cord

Powers LIFEPAK 20/20e defibrillator/monitors and recharges the internal battery.

- 11140-000015** (U.S.)
- 11140-000019** (Europe)
- 11140-000023** (South Africa)
- 11140-000045** (Unterminated)
- 11140-000017** (Australia)
- 11140-000022** (U.K.)
- 11140-000045** (Switzerland)



Rechargeable Internal Battery

For use as AC power backup.

- 11141-000112** (LIFEPAK 20e Lithium-ion)
- 11141-000068** (LIFEPAK 20 NiMH)



Lithium-ion Battery for CodeManagement Module™

Rechargeable battery for the CodeManagement Module; for use as AC power backup.

- 11141-000162** (Lithium-ion)

ECG Monitoring Accessories



5-Wire ECG Cable
11110-000066



3-Wire ECG Cable
11110-000029



LIFE-PATCH® ECG Electrodes
Adult, pregelled.
11100-000001 (3/pack)
11100-000002 (4/pack)

Cleartrace® ECG Electrodes
Adult, clear tape.
Available in the EU.
1700-003 (3/pack)



50mm ECG Printer Paper
11240-000013 (3 rolls/box)



CodeManagement Module for use with LIFEPAK 20/20e Defibrillator/Monitor - Wireless

Connects to the LIFEPAK 20 or 20e defibrillator to add wireless connectivity
11150-000018

CodeManagement Module for use with LIFEPAK 20/20e Defibrillator/Monitor - Wireless & Capnography

Connects to the LIFEPAK 20 or 20e defibrillator to add wireless connectivity and capnography monitoring
11150-000019

Therapy Delivery Accessories

Hard Paddles and Electrode Gel



Standard Adult Detachable Hard Paddles

Adult paddle plates are removable to reveal pediatric paddles.

11130-000037 (English)	11130-000052 (Polish)
11130-000038 (French)	11130-000053 (Brazilian)
11130-000045 (German/Dutch)	11130-000054 (Japanese)
11130-000046 (Spanish)	11130-000055 (Chinese)
11130-000047 (Italian)	11130-000057 (Hungarian)
11130-000048 (Swedish)	11130-000058 (Czech)
11130-000049 (Danish)	11130-000059 (Russian)
11130-000050 (Finnish)	11130-000060 (Korean)
11130-000051 (Norwegian)	21330-001024 (Replacement Plate)



SIGNAGEL® Electrode Gel

For use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel.

21300-005847 (8.5 oz/tube)



Internal Paddle Handles with Discharge Control

For use with LIFEPAK 15, 12 or 20/20e defibrillator/monitors.

11131-000001 (1 pair)



Internal Paddles (1 pair)

Requires Internal Paddle Handles

1" size

11131-000010 (6.25" shaft)

1.5" size

11131-000011 (6" shaft)

11131-000021 (9" shaft)

11131-000024 (14" shaft)

2" size

11131-000012 (5.75" shaft)

11131-000022 (8.75" shaft)

2.5" size

11131-000013 (5.5" shaft)

11131-000019 (8.5" shaft)

3.5" size

11131-000014 (5" shaft)

11131-000023 (8" shaft)

Therapy Delivery Accessories

EDGE System™ Electrodes for Pacing/Defibrillation/ECG with QUIK-COMBO® Connector

18 month minimum shelf life unless specified



EDGE System Electrodes with QUIK-COMBO Connector

11996-000091 (24" leadwire length)



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

11996-000090 (24" leadwire length)



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators; 12 month minimum shelf life at time of shipment.

11996-000093 (24" leadwire length)



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System

11996-000017 (42" leadwire length)



QUIK-COMBO Therapy Cable

11113-000040



FAST-PATCH® Adapter Cable

For use with the QUIK-COMBO therapy cable.

11110-000052

Other EDGE System™ Electrodes for Pacing/Defibrillation/ECG



EDGE System Electrodes with FAST-PATCH Connector

Must be used with QUIK-COMBO to FAST-PATCH Adapter Cable (11110-000052).

11996-000092

End-Tidal CO₂ (EtCO₂) Monitoring Accessories

Third Party Accessories



Microcap® Handheld Capnograph

Includes rechargeable battery pack, one Adult/Pediatric FilterLine® SET, one Adult Nasal FilterLine, one Adult Airway Adapter, one Infant/Neonatal Airway Adapter and operating instructions.

11996-000256

Microcap Plus Handheld Capnograph and Pulse Oximeter

Includes rechargeable battery pack, one Adult/Pediatric FilterLine SET, one Adult Nasal FilterLine, one ASP-3 Oxisensor II Adult Sensor Assortment Pack, one Adult Airway Adapter, one Infant/Neonatal Airway Adapter and operating instructions.

11996-000258

End-Tidal CO₂ (EtCO₂) Monitoring Accessories

Oridion® Filterlines for Intubated Patients

Single patient use



FilterLine® SET

Key applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

11996-000081 (25/pack, 200 cm)

11996-000164 (25/pack, 400 cm)



FilterLine H SET

Key applications: Critical Care, Humidified Environments.

Adult/Pediatric

11996-000080 (25/pack, 200 cm)

Adult/Neonatal

11996-000001 (25/pack, 200 cm)

Oridion Non-Intubated Filterlines

Single patient use



Smart CapnoLine® Plus

Key applications: Procedural sedation, Lower GI procedures, MAC, EMS, ED, Rapid Response Teams.

Adult with O₂

11996-000163 (25/pack, 200 cm)

11996-000167 (100/pack, 200 cm)

11996-000165 (25/pack, 400 cm)

Adult without O₂

11996-000162 (25/pack, 200 cm)

11996-000166 (100/pack, 200 cm)



Smart CapnoLine

Key applications: Procedural sedation, EMS, ED, Rapid Response Teams.

Pediatric with O₂

11996-000128 (25/pack, 200 cm)

Pediatric without O₂

11996-000120 (25/pack, 200 cm)

SpO₂ Monitoring Accessories

Masimo LNC® Reusable Patient Cables



LNC Patient Cable

- 11171-000024** (4 ft)
- 11171-000016** (10 ft)
- 11171-000025** (14 ft)
- 11171-000027** (4 ft extension)

Masimo LNC Compatible SpO₂ Only Sensors



LNC Adhesive Sensor

Single Patient Use. (20/box)

- 11171-000019** (Ad)
- 11171-000019** (Ped)
- 11171-000029** (Neo/Pt)
- 11171-000028** (Neo/Ad)
- 11171-000031** (Inf)



LNC Reusable Sensor

- 11171-000017** (Ad)
- 11171-000018** (Ped)



LNC Reusable Soft Sensor

- 11171-000052** (Ad)

Masimo Reusable Direct Connect Sensors and Cables



Direct Connect Reusable Soft Sensor

- 11171-000051** (Ad)

SpO₂ Monitoring Accessories

Masimo® LNOP Patient Cables



LNOP Patient Cable

11171-000006 (4 ft)

11171-000008 (8 ft)

11171-000009 (12 ft)

Masimo LNOP Patient Cable Compatible SpO₂ Only Sensors



LNOP Reusable Sensor

11171-000007 (Ad)

11171-000010 (Ped)



LNOP Adhesive Sensor (20/box)

11171-000011 (Ad)

11171-000012 (Ped)



LNOP Adhesive Sensor (20/box)

11171-000014 (Neo/Pt)

11171-000013 (Neo)

11171-000036 (Inf)

Additional Accessories and Adapter Cables



Reusable Ambient Light Shield

11171-000054 (5/pack)



Disposable Ambient Light Shield

11171-000055 (10/pack)

SpO₂ Monitoring Accessories

Masimo to Nellcor Adapter



MNC-1 Adapter Cable

Allows specially configured LIFEPAK 20e defibrillator/monitors with Masimo SpO₂ to connect to Nellcor sensors

11996-000183 (10 ft)

11996-000198 (4 ft)

Masimo to Nellcor Adapter Compatible SpO₂ Sensors



DURASENSOR Reusable Clip

11996-000060 (Ad)

DURA-Y Multisite Reusable Sensor

11996-000106 (≥ 1 kg)



Oxiband Reusable Sensor

Includes 50 disposable adhesive sensors

11996-000061 (Ad/Neo)

11996-000062 (Ped/Inf)

Disposable Adhesive Bandage Wrap (100/pack)

Not Available in Canada.

11996-000048 (Ad/Neo)

11996-000049 (Ped/Inf)



Oxisensor II Adhesive Sensors (24/box)

11996-000113 (Ad, 18 in.)

11996-000114 (Ad, 36 in.)

11996-000116 (Ped, 18 in.)



Oxisensor II Adhesive Sensors (24/box)

11996-000115 (Inf, 18 in.)

11996-000117 (Neo/Ad 18 in.)

Carrying Cases and Mounting Options



Basic Carrying Case/ Accessories Organizer

Includes set of right and left pouches for additional accessories.

11260-000018



Carrying Case for LIFEPAK 20/20e Defibrillator with Module

Includes a set of right and left pouches for additional storage of accessories. For use with a LIFEPAK 20/20e device with CodeManagement Module.

11260-000045



Top Pouch

For additional storage: insert in place of standard paddles.

11260-000043

Shoulder Strap for Basic Carrying Case

Includes mounting brackets and hardware.

11260-000041



Docking Station

Attaches to crash cart to lock LIFEPAK 20/20e defibrillator/monitor safely in place.

21330-000996

Right Side Accessory Pouch

Attaches to LIFEPAK 20/20e defibrillator/monitors and holds items such as QUIK-COMBO therapy cable, electrodes and SIGNAGEL electrode gel.

11260-000016

Training Tools and Testers



Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias and bradycardias, as well as ST segment and T wave abnormalities.

11996-000311 (12-Lead)

11996-000310 (3-Lead)



Reusable Training Electrodes

For use with QUIK-COMBO patient simulators. Allows users to practice defibrillation electrode placement. Includes instructions for use and reusable training electrodes.

11250-000052 (5/pair)

11101-000003 (Replacement 5/pair)



Defibrillator Checker

Tests integrity of energy delivery through the standard hard paddles. Neon light indicates energy has been delivered.

11998-000060



QUIK-COMBO Test Plug

Test plug connects to QUIK-COMBO therapy cable to test cable integrity.

11113-000002

Literature, Communication and Third Party Accessories



Inservice DVD

3314554-000 (English - DVD)

Service Manual

26500-003543 (English - CD - LIFEPAK 20/20e with CodeManagement Module)

Operating Instructions

26500-003224 (Intl English)
26500-003223 (English)
26500-003235 (Brazilian)
26500-003233 (Chinese)
26500-003244 (Czech)
26500-003227 (Danish)
26500-003237 (Dutch)
26500-003228 (Finnish)
26500-003241 (French)
26500-003242 (German)
26500-003226 (Hungarian)
26500-003234 (Italian)
26500-003232 (Korean)
26500-003229 (Norwegian)
26500-003230 (Polish)
26500-003239 (Portuguese)
26500-003248 (Russian)
26500-003240 (Spanish)
26500-003236 (Swedish)

Communication Accessories



Serial Port Cable

For connecting LIFEPAK 20/20e defibrillator/monitor to PC.

11230-000018 (6 ft)

Transport Configuration Cable

For connecting LIFEPAK 20/20e to LIFEPAK 20/20e defibrillator/monitor.

11230-000019 (6 ft)



Physio-Control Headquarters
11811 Willows Road NE
Redmond, WA 98052 USA
Tel 425 867 4000
Fax 425 867 4121
www.physio-control.com

**Physio-Control Operations
Netherlands B.V.**
Keizersgracht 125-127,
1015 CJ Amsterdam, NL
Tel +31 (0)20 7070560
Fax +31 (0)20 3301194
www.physio-control.nl

**Physio-Control
UK Sales Ltd**
Concorde House
Trinity Park, Solihull
Birmingham
B37 7UQ
United Kingdom
Tel (0)44 1923 212 213
Fax (0)44 1923 241 004

**Physio-Control
Australia Pty Ltd**
Suite 4.01
15 Orion Road
Lane Cove
NSW 2066
Australia
Toll Free Tel 1800 987 982
Toll Free Fax 1800 890 892

For further information please contact your local Physio-Control representative or visit our website at www.physio-control.com