



# LIFEPAK<sup>®</sup> 20

DEFIBRILLATOR/MONITOR SERIES

Genuine Accessories from Physio-Control

# LIFEPAK 20 MONITOR/DEFIBRILLATOR SERIES





## Genuine Accessories from Physio-Control

Ensure the safety of your staff and patients. Only Physio-Control accessories and disposables have been thoroughly tested with LIFEPAK products. Stick with the name you trust and the company that stands behind it. To order, contact your Physio-Control representative.

## Power Options

---



### AC Power Cord

Powers LIFEPAK 20/20e defibrillator/monitors and recharges the internal battery.

- 11140-000015** (U.S.)
- 11140-000019** (Europe)
- 11140-000023** (South Africa)
- 11140-000045** (Unterminated)
- 11140-000017** (Australia)
- 11140-000022** (U.K.)
- 11140-000045** (Switzerland)



### Rechargeable Internal Battery

For use as AC power backup.

- 11141-000112** (LIFEPAK 20e Lithium-ion)
- 11141-000068** (LIFEPAK 20 NiMH)



### Lithium-ion Battery for CodeManagement Module™

Rechargeable battery for the CodeManagement Module; for use as AC power backup.

- 11141-000162** (Lithium-ion)

# ECG Monitoring Accessories



**5-Wire ECG Cable**  
11110-000066



**3-Wire ECG Cable**  
11110-000029



**LIFE-PATCH® ECG Electrodes**  
Adult, pregelled.  
11100-000001 (3/pack)  
11100-000002 (4/pack)

**Cleartrace® ECG Electrodes**  
Adult, clear tape.  
Available in the EU.  
1700-003 (3/pack)



**50mm x 22mm ECG Printer Paper**  
11240-000013 (3 rolls/box)



**CodeManagement Module for use with LIFEPAK 20/20e Defibrillator/Monitor - Wireless**

Connects to the LIFEPAK 20 or 20e defibrillator to add wireless connectivity  
11150-000018

**CodeManagement Module for use with LIFEPAK 20/20e Defibrillator/Monitor - Wireless & Capnography**

Connects to the LIFEPAK 20 or 20e defibrillator to add wireless connectivity and capnography monitoring  
11150-000019

# Therapy Delivery Accessories

## Hard Paddles and Electrode Gel



### Standard Adult Detachable Hard Paddles

Adult paddle plates are removable to reveal pediatric paddles.

<b>11130-000037</b> (English)	<b>11130-000052</b> (Polish)
<b>11130-000038</b> (French)	<b>11130-000053</b> (Brazilian)
<b>11130-000045</b> (German/Dutch)	<b>11130-000054</b> (Japanese)
<b>11130-000046</b> (Spanish)	<b>11130-000055</b> (Chinese)
<b>11130-000047</b> (Italian)	<b>11130-000057</b> (Hungarian)
<b>11130-000048</b> (Swedish)	<b>11130-000058</b> (Czech)
<b>11130-000049</b> (Danish)	<b>11130-000059</b> (Russian)
<b>11130-000050</b> (Finnish)	<b>11130-000060</b> (Korean)
<b>11130-000051</b> (Norwegian)	<b>21330-001024</b> (Replacement Plate)



### SIGNAGEL® Electrode Gel

For use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel.

**21300-005847** (8.5 oz/tube)



### Internal Paddle Handles with Discharge Control

For use with LIFEPAK 15, 12 or 20/20e defibrillator/monitors.

**11131-000001** (1 pair)



### Internal Paddles (1 pair)

Requires Internal Paddle Handles

#### 1" size

**11131-000010** (6.25" shaft)

#### 1.5" size

**11131-000011** (6" shaft)

**11131-000021** (9" shaft)

**11131-000024** (14" shaft)

#### 2" size

**11131-000012** (5.75" shaft)

**11131-000022** (8.75" shaft)

#### 2.5" size

**11131-000013** (5.5" shaft)

**11131-000019** (8.5" shaft)

#### 3.5" size

**11131-000014** (5" shaft)

**11131-000023** (8" shaft)

# Therapy Delivery Accessories

## EDGE System™ Electrodes for Pacing/Defibrillation/ECG with QUIK-COMBO® Connector

18 month minimum shelf life unless specified



EDGE System Electrodes with QUIK-COMBO Connector

11996-000091 (24" leadwire length)



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

11996-000090 (24" leadwire length)



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators;  
12 month minimum shelf life at time of shipment.

11996-000093 (24" leadwire length)



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System

11996-000017 (42" leadwire length)



QUIK-COMBO Therapy Cable

11110-000040



FAST-PATCH® Adapter Cable

For use with the QUIK-COMBO therapy cable.  
11110-000052

## Other EDGE System™ Electrodes for Pacing/Defibrillation/ECG



EDGE System Electrodes with FAST-PATCH Connector

Must be used with QUIK-COMBO to FAST-PATCH Adapter Cable (11110-000052).

11996-000092

# End-Tidal CO<sub>2</sub> (EtCO<sub>2</sub>) Monitoring Accessories

---

## Oridion<sup>®</sup> Filterlines for Intubated Patients

Single patient use



### FilterLine<sup>®</sup> SET

Key applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

**11996-000081** (25/pack, 200 cm)

**11996-000164** (25/pack, 400 cm)



### FilterLine H SET

Key applications: Critical Care, Humidified Environments.

Adult/Pediatric

**11996-000080** (25/pack, 200 cm)

Adult/Neonatal

**11996-000001** (25/pack, 200 cm)

---

## Oridion Non-Intubated Filterlines

Single patient use



### Smart CapnoLine<sup>®</sup> Plus

Key applications: Procedural sedation, Lower GI procedures, MAC, EMS, ED, Rapid Response Teams.

Adult with O<sub>2</sub>

**11996-000163** (25/pack, 200 cm)

**11996-000167** (100/pack, 200 cm)

**11996-000165** (25/pack, 400 cm)

Adult without O<sub>2</sub>

**11996-000162** (25/pack, 200 cm)

**11996-000166** (100/pack, 200 cm)



### Smart CapnoLine

Key applications: Procedural sedation, EMS, ED, Rapid Response Teams.

Pediatric with O<sub>2</sub>

**11996-000128** (25/pack, 200 cm)

Pediatric without O<sub>2</sub>

**11996-000120** (25/pack, 200 cm)



# SpO<sub>2</sub> Monitoring Accessories

---

## Masimo LNC® Reusable Patient Cables



### LNC Patient Cable

**11171-000024** (4 ft)  
**11171-000016** (10 ft)  
**11171-000025** (14 ft)  
**11171-000027** (4 ft extension)

## Masimo LNC Compatible SpO<sub>2</sub> Only Sensors



### LNC Adhesive Sensor

Single Patient Use. (20/box)

**11171-000019** (Ad)  
**11171-000020** (Ped)  
**11171-000029** (Neo/Pt)  
**11171-000028** (Neo/Ad)  
**11171-000031** (Inf)



### LNC Reusable Sensor

**11171-000017** (Ad)  
**11171-000018** (Ped)



### LNC Reusable Soft Sensor

**11171-000052** (Ad)

## Masimo Reusable Direct Connect Sensors and Cables



### Direct Conect Reusable Soft Sensor

**11171-000051** (Ad)

## SpO<sub>2</sub> Monitoring Accessories

---

### Masimo® LNOP Patient Cables



LNOP Patient Cable

11171-000006 (4 ft)

11171-000008 (8 ft)

11171-000009 (12 ft)

---

### Masimo LNOP Patient Cable Compatible SpO<sub>2</sub> Only Sensors



LNOP Reusable Sensor

11171-000007 (Ad)

11171-000010 (Ped)



LNOP Adhesive Sensor (20/box)

11171-000011 (Ad)

11171-000012 (Ped)



LNOP Adhesive Sensor (20/box)

11171-000014 (Neo/Pt)

11171-000013 (Neo)

11171-000036 (Inf)

---

### Additional Accessories and Adapter Cables



Reusable Ambient Light Shield

11171-000054 (5/pack)



Disposable Ambient Light Shield

11171-000055 (10/pack)

# SpO<sub>2</sub> Monitoring Accessories

## Masimo to Nellcor Adapter



### MNC-1 Adapter Cable

Allows specially configured LIFEPAK 20e defibrillator/monitors with Masimo SpO<sub>2</sub> to connect to Nellcor sensors

**11996-000183** (10 ft)

**11996-000198** (4 ft)

## Masimo to Nellcor Adapter Compatible SpO<sub>2</sub> Sensors



### DURASENSOR Reusable Clip

**11996-000060** (Ad)

### DURA-Y Multisite Reusable Sensor

**11996-000106** (≥ 1 kg)



### Oxiband Reusable Sensor

Includes 50 disposable adhesive sensors

**11996-000061** (Ad/Neo)

**11996-000062** (Ped/Inf)

### Disposable Adhesive Bandage Wrap (100/pack)

Not Available in Canada.

**11996-000048** (Ad/Neo)

**11996-000049** (Ped/Inf)



### Oxisensor II Adhesive Sensors (24/box)

**11996-000113** (Ad, 18 in.)

**11996-000114** (Ad, 36 in.)

**11996-000116** (Ped, 18 in.)



### Oxisensor II Adhesive Sensors (24/box)

**11996-000115** (Inf, 18 in.)

**11996-000117** (Neo/Ad 18 in.)

## Carrying Cases and Mounting Options

---



### Basic Carrying Case/ Accessories Organizer

Includes set of right and left pouches for additional accessories.  
**11260-000018**



### Carrying Case for LIFEPAK 20/20e Defibrillator with Module

Includes a set of right and left pouches for additional storage of accessories. For use with a LIFEPAK 20/20e device with CodeManagement Module.  
**11260-000045**



### Top Pouch

For additional storage: insert in place of standard paddles.  
**11260-000043**

### Shoulder Strap for Basic Carrying Case

Includes mounting brackets and hardware.  
**11260-000041**



### Docking Station

Attaches to crash cart to lock LIFEPAK 20/20e defibrillator/monitor safely in place.  
**21330-000996**

### Right Side Accessory Pouch

Attaches to LIFEPAK 20/20e defibrillator/monitors and holds items such as QUIK-COMBO therapy cable, electrodes and SIGNAGEL electrode gel.  
**11260-000016**

# Training Tools and Testers



## Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias and bradycardias, as well as ST segment and T wave abnormalities.

**11996-000311** (12-Lead)

**11996-000310** (3-Lead)



## Defibrillator Checker

Tests integrity of energy delivery through the standard hard paddles. Neon light indicates energy has been delivered.

**11998-000060**



## QUIK-COMBO Test Plug

Test plug connects to QUIK-COMBO therapy cable to test cable integrity.

**11113-000002**

# Literature, Communication and Third Party Accessories

---



## Inservice DVD

**26500-003545** (English - DVD)

## Service Manual

**26500-003543** (English - CD - LIFEPAK 20/20e with CodeManagement Module)

## Operating Instructions

**26500-003224** (Intl English)  
**26500-003223** (English)  
**26500-003235** (Brazilian)  
**26500-003233** (Chinese)  
**26500-003244** (Czech)  
**26500-003227** (Danish)  
**26500-003237** (Dutch)  
**26500-003228** (Finnish)  
**26500-003241** (French)  
**26500-003242** (German)  
**26500-003226** (Hungarian)  
**26500-003234** (Italian)  
**26500-003232** (Korean)  
**26500-003229** (Norwegian)  
**26500-003230** (Polish)  
**26500-003239** (Portuguese)  
**26500-003248** (Russian)  
**26500-003240** (Spanish)  
**26500-003236** (Swedish)

---

## Communication Accessories



## Serial Port Cable

For connecting LIFEPAK 20/20e defibrillator/monitor to PC.

**11230-000018** (6 ft)

## Transport Configuration Cable

For connecting LIFEPAK 20/20e to LIFEPAK 20/20e defibrillator/monitor.

**11230-000019** (6 ft)



For yderligere information kontakt veligst deres lokale Physio-Control representant eller besøg vores hjemmeside [www.physio-control.com](http://www.physio-control.com)



**Physio-Control Headquarters**  
11811 Willows Road NE  
Redmond, WA 98052 USA  
Tel 425 867 4000  
Fax 425 867 4121  
[www.physio-control.com](http://www.physio-control.com)

**Physio-Control Operations  
Netherlands B.V.**  
HQ Europe, Middle East & Africa  
UP Office Building  
Piet Heinkade 55,  
1019 GM Amsterdam NL  
Tel +31 (0)20 7070560  
Fax +31 (0)20 3391194  
[www.physio-control.com](http://www.physio-control.com)

**Physio-Control Denmark Sales ApS**  
[info.nordic@physio-control.com](mailto:info.nordic@physio-control.com)  
Tel 69 60 70 58  
[www.physio-control.dk](http://www.physio-control.dk)