

McGRATH[®]
MAC EMS

Video Laryngoscope

ACCESSORIES

Designed & manufactured
by Aircraft Medical



Distributed by **Physio-Control**

McGRATH[®] MAC EMS

Video Laryngoscope





**Thank you for choosing Physio-Control
as your partner in helping save lives and
improve patient care.**

The McGRATH® MAC EMS video laryngoscope is shipped with one battery and the Instructions for Use. Disposable blades are shipped separately. In this catalog you will find an overview of blades and additional accessories available for the McGRATH® MAC EMS. If you have any questions about accessories or disposables, please contact your Physio-Control representative.

Disposable Laryngoscope Blades

McGRATH® MAC Disposable Laryngoscope Blades are high performance, low-cost, single-use sterile blades. They are used for routine as well as many difficult airways, from pediatric to adult patients. They have a Macintosh-like curvature and are therefore easy to use with limited training. McGRATH® MAC Blades are made of robust optical polymer and blade lenses are coated with hydrophilic optical surface treatment. They are packed individually and shipped in cartons of 10 or 50.



McGRATH® MAC 2 Disposable Laryngoscope Blade

Designed for pediatric patients older than eight weeks or at least 4.5 kg. Packed individually and shipped in cartons of 10 or 50.

Carton of 50 11996-000395

Carton of 10 11996-000414



McGRATH® MAC 3 Disposable Laryngoscope Blade

Designed for adult patients. Packed individually and shipped in cartons of 10 or 50.

Carton of 50 11996-000396

Carton of 10 11996-000415



McGRATH® MAC 4 Disposable Laryngoscope Blade


Designed for large adult patients. Packed individually and shipped in cartons of 10 or 50.

Carton of 50 11996-000397

Carton of 10 11996-000416

McGRATH® X blade™

The McGRATH® X blade™ is a sterile disposable single-use blade that extends the range of proprietary disposable laryngoscope blades for the McGRATH® MAC EMS system. The X blade™ contains a sweeping acute curvature for extreme airways, complementing the existing range of McGRATH® MAC Disposable Laryngoscope Blades which are utilized in routine to difficult airways.



McGRATH® X3 Disposable Laryngoscope Blade
 Designed for adult patients with extreme anterior airways.
 Packed individually and shipped in cartons of 10.
11996-000398

Replacement Batteries



The McGRATH® MAC EMS utilizes 3.6V Lithium power cells, which provide up to 250 minutes of usage. A representation of minutes remaining is displayed on screen. The on/off power switch is integrated into the battery pack, so when a new battery is fitted, the on/off function is renewed, keeping the product young.

McGRATH® 3.6V EMS Battery

The McGRATH® MAC EMS is shipped with one 3.6V EMS Battery. Replacement batteries are shipped individually.
11996-000394

McGRATH® MAC EMS

Video Laryngoscope





All accessories are latex free.

For further information please contact your local Physio-Control representative or visit our website at www.physio-control.com.



Physio-Control Headquarters
11811 Willows Road NE
Redmond, WA 98052 USA
Tel 425 867 4000
Fax 425 867 4121
www.physio-control.com

**Physio-Control Operations
Netherlands B.V.**
HQ Europe, Middle East & Africa
UP Office Building
Piet Heinkade 55,
1019 GM Amsterdam NL
Tel +31 (0)20 7070560
Fax +31 (0)20 3391194
www.physio-control.com

**Physio-Control
UK Sales Ltd**
Concorde House
Trinity Park, Solihull
Birmingham
B37 7UQ
United Kingdom
Tel +44 (0)808 258 0094
Fax +44 (0)808 189 0337

**Physio-Control
Australia Pty Ltd**
Suite 4.01
15 Orion Road
Lane Cove
NSW 2066
Australia
Toll Free Tel 1800 987 982
Toll Free Fax 1800 890 892



Aircraft Medical Limited, 7E Cross Way, Dalgety Bay, United Kingdom KY11 9JE

