



LIFEPAK[®] 15

MONITOR/DEFIBRILLATOR

Genuine Accessories from Physio-Control

LIFEPAK[®] 15 MONITOR/DEFIBRILLATOR



Genuine Accessories from Physio-Control.

Ensure the safety of your staff and patients. Only Physio-Control accessories and disposables have been thoroughly tested with LIFEPAK products. Stick with the name you trust and the company that stands behind it. To order contact your Physio-Control representative.

Power Options



REDI-CHARGE Base

AC power cord and adapter tray not included.

11141-000115

REDI-CHARGE LIFEPAK 15 Adapter Tray

11140-000052

AC Power Cords

11140-000015 (U.S.)

11140-000019 (Europe)

11140-000023 (South Africa)

11140-000046 (Italian)

11140-000017 (Australia)

11140-000021 (U.K.)

11140-000045 (Switzerland)



Mobile Battery Charger

Includes AC and DC power cords, mounting bracket and operating instructions.

11577-000011 (U.S.)

11577-000015 (United Kingdom)

11577-000012 (Europe)

11577-000017 (Australia)

11577-000014 (Switzerland)



Rechargeable Lithium-ion Battery

5.7 ah, 10.8 volt.

21330-001176



AC Power Adapter

Includes Right Angle Cable (AC power cord not included)

11140-000072

DC Power Adapter

Includes DC cable and Right Angle Cable

11140-000074

AC Power Cord

11140-000015 (U.S.)

11140-000019 (Europe)

11140-000023 (South Africa)

11140-000046 (Italian)

11140-000017 (Australia)

11140-000021 (U.K.)

11140-000045 (Switzerland)

Power Attachment Kit

11577-000019



Extension Cable

For AC/DC Power Adapter.

11140-000080



Replacement Right Angle Power Cable

For AC/DC Power Adapter.

11140-000081




Replacement DC Input Cable

For DC Power Adapter.

11140-000084

ECG Monitoring Accessories



12-Lead ECG Cable Trunk Cable with 4-Wire Limb Leads

11111-000018 (5ft, AHA)
11111-000019 (5ft, IEC)
11111-000020 (8ft, AHA)
11111-000021 (8ft, IEC)



12-Lead ECG Cable 6-Wire Precordial Attachment

11111-000022 (AHA)
11111-000023 (IEC)



4-wire Cable Comb

21300-008054 (10/pack)

6-wire Cable Comb

21330-008055 (10/pack)



3-Wire ECG Cable

11110-000029 (AHA)
11110-000030 (IEC)



5-Wire ECG Cable

4-wire limb plus 1 chest lead, labeled "V1" on the monitor reports.
11110-000066 (AHA)
11110-000067 (IEC)



Chart Recorder Paper

100mm x 22m
11240-000016 (2 rolls per box)



LIFE•PATCH® ECG Electrodes

Adult, pregelled
(Not available in the European Union)

11100-000001 (3/pack)
11100-000002 (4/pack)

Cleartrace® ECG Electrodes

Adult, clear tape
(Available in the European Union)

11100-000009 (3/pack)

Therapy Delivery Accessories

Hard Paddles and Electrode Gel



Standard Hard Paddles (1 pair)

- 11130-000061 (English)
- 11130-000062 (French)
- 11130-000063 (German)
- 11130-000064 (Spanish)
- 11130-000065 (Italian)
- 11130-000066 (Swedish)
- 11130-000067 (Danish)
- 11130-000068 (Portuguese)
- 11130-000069 (Norwegian)
- 11130-000070 (Dutch)
- 11130-000071 (Finnish)
- 11130-000072 (Polish)
- 11130-000073 (Japanese)



SIGNAGEL® Electrode Gel

For use with hard paddles. Recommended for ECG, defibrillation, biofeedback and EMG.
21300-005847 (8.5 oz)



Pediatric Paddle Adapter

11133-000007 (1 adapter)



Internal Paddles (Requires Internal Paddle Handles and Internal Paddles Adapter Cable)

- | | |
|--|---|
| 1 in. size | 2.5 in. size |
| 11131-000010 (1 pair, 6.25 in. shaft) | 11131-000013 (1 pair, 5.5 in. shaft) |
| 1.5 in. size | 11131-000019 (1 pair, 8.5 in. shaft) |
| 11131-000011 (1 pair, 6 in. shaft) | 3.5 in. size |
| 11131-000021 (1 pair, 9 in. shaft) | 11131-000014 (1 pair, 5 in. shaft) |
| 11131-000024 (1 pair, 14 in. shaft) | 11131-000023 (1 pair, 8 in. shaft) |
| 2 in. size | |
| 11131-000012 (1 pair, 5.75 in. shaft) | |
| 11131-000022 (1 pair, 8.75 in. shaft) | |



Internal Paddle Handles with Discharge Control

(For use with the Internal Paddles Adapter Cable)

11131-000001 (1 pair)



Internal Paddles Adapter Cable

(For use with Internal Paddle Handles)

11998-000326

Therapy Delivery Accessories

EDGE System™ Electrodes for Pacing/Defibrillation/ECG with QUIK-COMBO® Connector

18-month minimum shelf life remaining at time of shipment from Physio-Control except where noted.



EDGE System Electrodes with QUIK-COMBO Connector

24" leadwire length

11996-000091



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector

24" leadwire length

11996-000090



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK™ Preconnect System

42" leadwire length

11996-000017



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators;
12 month minimum shelf life at time of shipment
24" leadwire length.

11996-000093



QUIK-COMBO Therapy Cable

With convenient TRUE-LOCK™ Cable Connector.

11113-000004 (8 ft)

NIBP Monitoring Accessories

NIBP Hoses



NIBP Tubing

21300-008147 (9ft)
21300-008146 (12ft)



NIBP Tubing, Coiled

21300-008148 (2-9ft)

NIBP Cuffs



Reusable Cuff

X-Large Adult

35 - 44 cm
11160-000019

Pediatric

13 - 20 cm
11160-000013

Large Adult

32 - 42 cm
11160-000017

Infant

8 - 14 cm
11160-000011

Adult

26 - 35 cm
11160-000015



Single Patient Use Cuff

X-Large Adult

35 - 44 cm
11160-000020

Pediatric

13 - 20 cm
11160-000014

Large Adult

32 - 42 cm
11160-000018

Infant

8 - 14 cm
11160-000012

Adult

26 - 35 cm
11160-000016

Pulse Oximetry Monitoring Accessories

Masimo SET® RC Patient Cables



RC Patient Cable

For use with M-LNCS and Rainbow Patient Sensors

11171-000037 (4ft)

11171-000038 (12ft)

Masimo SET RC Patient Cable Compatible SpO₂ Sensors



M-LNCS Reusable Sensor

11171-000046 (Ad)

11171-000047 (Ped)



M-LNCS Adhesive Sensors (20/box)

11171-000039 (Ad)

11171-000040 (Ped)



M-LNCS Adhesive Sensors (20/box)

11171-000043 (Neo/Pt)

11171-000042 (Neo/Ad)

11171-000041 (Inf)

Masimo SET RC Patient Cable Compatible Rainbow® SpO₂, SpCO, SpMet Sensors



Rainbow Reusable Sensor

11171-000049 (Ad)

11171-000050 (Ped)



Rainbow Adhesive Sensor (10/box)

11996-000339 (Ad)

11996-000340 (Ped)



Rainbow Adhesive Sensor (10/box)

11996-000342 (Inf)

11996-000341 (Neo/Ad)

Masimo SET LNC Patient Cables



Red LNC Patient Cable

For Use with LNCS Patient Sensors

11996-000323 (4ft)

11996-000324 (10ft)

11996-000325 (14ft)

Masimo SET LNC Patient Cable Compatible SpO₂ Sensors



LNCS® Reusable Sensor

11171-000017 (Ad)

11171-000018 (Ped)



LNCS Reusable Soft Sensor

11171-000052 (Ad)



LNCS Adhesive Sensor (20/box)

11171-000019 (Ad)

11171-000020 (Ped)



LNCS Adhesive Sensor (20/box)

11171-000029 (Neo/Pt)

11171-000028 (Neo/Ad)

11171-000031 (Inf)

Pulse Oximetry Monitoring Accessories

Direct Connect SpO₂ Only Patient Sensors



Adult Reusable Direct Connect Sensor
11996-000331 (3ft)
11996-000332 (12ft)



Pediatric Reusable Direct Connect Sensor
11996-000333 (3ft)
11996-000334 (12ft)



Adult Reusable Soft Direct Connect Sensor
11171-000053 (8ft)

Direct Connect Rainbow SpO₂, SpCO, SpMet Patient Sensors



Adult Rainbow Direct Connect Reusable Sensor
11996-000335 (3ft)
11171-000032 (8ft)
11996-000336 (12ft)



Pediatric Rainbow Direct Connect Reusable Sensor
11996-000337 (3ft)
11171-000033 (8ft)
11996-000338 (12ft)

Additional Masimo Accessories



Reusable Ambient Light Shield
11171-000054 (5/bag)



Disposable Ambient Light Shield
11171-000055 (10/bag)

Pulse Oximetry Monitoring Accessories

Masimo to Nellcor Adapter



Red MNC Cable

Connects LIFEPAK 15 to Nellcor patient sensor

11996-000365 (4ft)

11996-000366 (10ft)

Masimo to Nellcor Adapter Compatible SpO₂ Sensors



DURASENSOR Reusable Clip

11996-000060 (Ad)

DURA-Y Multisite Reusable Sensor

11996-000106 (> 1 kg)



Oxiband Reusable Sensor

Includes 50 disposable adhesive sensors.

11996-000061 (Ad/Neo)

11996-000062 (Ped/Inf)

Disposable Adhesive Bandage Wrap (100/pk)

Not Available in Canada.

11996-000048 (Ad/Neo)

11996-000049 (Ped/Inf)

Temperature Monitoring

Temperature Cables



Temperature Adapter Cable

11140-000079 (10ft)

11140-000078 (5ft)

Temperature Sensors

Single patient use



Esophageal-Rectal Sensor

11996-000360 (9FR, 20/box)



Skin Probe Sensor

11996-000359 (20/box)



Foley Catheter Sensor

11996-000361 (14FR, 10/pk)

11996-000362 (16FR, 10/pk)

11996-000363 (18FR, 10/pk)

End-Tidal CO₂ (EtCO₂) Monitoring Accessories

Oridion® Filterlines for Intubated Patients

Single patient use



FilterLine® SET

Key Applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

11996-000081 (25/pk, 200 cm)

11996-000164 (25/pk, 400 cm)



FilterLine H SET

Key Applications: Critical care, Humidified Environments.

Adult/Pediatric

11996-000080 (25/pk, 200 cm)

Infant/Neonatal

11996-000001 (25/pk, 200 cm)

Oridion Non-Intubated Filterlines

Single patient use



Smart CapnoLine® Plus

Key Applications: Procedural Sedation, Upper GI Procedures, MAC, RMS, BD, Rapid Response Teams.

Adult with O₂

11996-000163 (25/pk, 200 cm)

11996-000167 (100/pk, 200 cm)

11996-000165 (25/pk, 400 cm)

Adult without O₂

11996-000162 (25/pk, 200 cm)

11996-000166 (100/pk, 200 cm)



Smart CapnoLine

Key Applications: Procedural Sedation, SMS, SD, Rapid Response Teams.

Pediatric with O₂

11996-000128 (25/pk, 200 cm)

Pediatric without O₂

11996-000120 (25/pk, 200 cm)

Cases & Mounting Options



Standard Carrying Case

Includes right pouch and left pouch.
11577-000002

Shoulder Strap

11577-000001



Bed Connector

For use in hospital only.
11996-000374

Top Pouch

Storage for sensors and electrodes; insert in place of standard paddles.
11220-000028

Back Pouch

Ideal for additional accessory storage.
11260-000039

Communication Accessories



LIFEPAK Monitor to PC Cable

For connecting LIFEPAK 12 or LIFEPAK 15 monitor/defibrillator to PC.

11230-000020 (Serial)
11996-000369 (USB)



Titan II Wireless Gateway

For transmitting data from LIFEPAK 12/15 to the LIFENET™ System or CODE-STAT™ data review software. Requires existing wireless network.

21996-000073 (Wireless Gateway)
21996-000093 (Wireless 3G Gateway - GSM)
21996-000092 (Wireless 3G + Voice Gateway - GSM)



3G Gateway

For transmitting data from LIFEPAK 12/15 to the LIFENET System or CODE-STAT data review software. Requires data plan.

21996-000086 (GSM)

Testers & Training Materials



Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias, and bradycardias, as well as ST segment and T wave abnormalities.

11996-000311 (12-Lead)

11996-000310 (3-Lead)



Defibrillator Checker

Tests integrity of energy delivery through standard hard paddles. Neon light indicates energy has been delivered.

11998-000060



Test Load

Use to perform therapy cable performance checks. Connects to the QUIK-COMBO therapy cable on the defibrillator.

21330-001368 (Brazilian, Hungarian, Portuguese, Polish, Romanian)

21330-001369 (Chinese, Czech, Japanese, Korean, Russian)

21330-001366 (Danish, Finnish, Norwegian, Swedish)

21330-001367 (Dutch, French Canadian, German, Intl English, Italian, Spanish)

21330-001365 (English)

Training Tools



Inservice Video:

21330-001357 (English - NTSC DVD)

Operating Instructions:

26500-003255 (English)

26500-003256 (Intl English)

26500-003330 (German)

26500-003331 (Italian)

26500-003332 (French)

26500-003333 (Dutch)

26500-003334 (Spanish)

26500-003335 (Iberian Portuguese)

26500-003336 (Brazilian Portuguese)

26500-003337 (Swedish)

26500-003338 (Danish)

26500-003339 (Finnish)

26500-003340 (Norwegian)

26500-003341 (Polish)

26500-003342 (Hungarian)

26500-003343 (Czech)

26500-003344 (Russian)

26500-003345 (Korean)

Physio-Control er nu en del af Stryker.

For yderligere information kontakt veligst deres lokale Physio-Control Repræsentant eller besøg vores hjemmeside www.physio-control.com

Physio-Control Headquarters

11811 Willows Road NE
Redmond, WA 98052 USA
Tel 425 867 4000
www.physio-control.com

Physio-Control Operations

Netherlands B.V.
Galjoenweg 68
6222 NV Maastricht
The Netherlands
Tel +31 (0)43 3620008
Fax +31 (0)43 3632001

Physio-Control Denmark Sales ApS

info.nordic@physio-control.com
Tel +45 69 60 70 58
www.physio-control.dk