

WALL & FREESTANDING CABINETS

Cabinets are a great way to securely store your AED in public.

Cabinets increase device visibility and allow for easy access in times of emergency. All cabinets are International Building Code certified and comfortably fit all models of Physio-Control AEDs. Alarm and strobe light features help ensure your AEDs are protected.

Color Choices



White



Stainless Steel



OSHA Orange



Surface Mount (Projects 7" from Wall)

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11998-000310	White	Carbon Steel	17.38" x 17.38" x 7"	Square	•	•	
11220-000076	Stainless Steel	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		
11220-000079	White	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		
11220-000082	OSHA Orange	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		
11220-000083	White	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•	•	
11220-000084	Stainless Steel	Carbon Steel	14" x 14" x 6.75"	Rolled	•	•	
11220-000086	OSHA Orange	Aluminum	17.38" x 17.38" x 7"	Square	•	•	
11220-000087	OSHA Orange	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•	•	



Semi-Recessed (Projects 3" from Wall)

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11998-000292	White	Carbon Steel	14" x 14" x 6.75"	Rolled	•		
11210-000026	White	Carbon Steel	14" x 14" x 6.75"	Rolled	•		•
11220-000077	Stainless Steel	Carbon Steel	14" x 14" x 6.75"	Rolled	•		



Recessed Mount (Projects 1.5" from Wall)

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11210-000027	White	Carbon Steel	14" x 14" x 6.75"	Square	•		•
11220-000078	Stainless Steel	Carbon Steel	14" x 14" x 6.75"	Square	•		
11220-000088	OSHA Orange	Aluminum	14" x 14" x 6.75"	Square	•	•	
11220-000089	OSHA Orange	Carbon Steel	14" x 14" x 6.75"	Square	•		
11220-000090	OSHA Orange	Carbon Steel	14" x 14" x 6.75"	Square	•	•	
11998-000293	White	Carbon Steel	14" x 14" x 6.75"	Square	•		



Free Standing

Cat. Number	Color	Material	Size	Edge	Alarm	Strobe	Fire Resistant
11210-000028	White	Carbon Steel	17.38" x 17.38" x 7"	Square	•		
11210-000029	Gray	Carbon Steel	17.38" x 17.38" x 7"	Rolled	•		

AED LOCATION WALL SIGNS

AED Location Wall Signs alert responders to the location of an AED.

AED Location Wall Signs come in various designs and models to match the environment while meeting local AED recommendations. Our new, lighter materials are less obtrusive and easier to hang while brighter designs improve visibility. Help decrease response times by signaling the location of your AED.



Tent AED Location Sign



Flat AED Location Sign



T-mount AED Location Sign

Cat. Number	Style	Material	Size Wide x Tall	Thickness	Weight	Depth from Wall	Width across Wall
11998-000332	Tent w/ Logo	White Styrene	11" x 8"	0.079"	4 oz	4" deep	11"
11998-000330	Flat w/ Logo	White Styrene	8.5" x 11"	0.02"	1 oz	Flush	8.5"
11998-000331	T-mount w/ Logo	White PVC	8" x 10"	0.125"	4 oz	8" deep	1" wide lip
11998-000333	T-mount w/o Logo	White PV	8" x 10"	0.125"	4 oz	8" deep	1" wide lip

For further information, please contact Physio-Control at 800.442.1142 (U.S.), 800.895.5896 (Canada) or visit our website at www.physio-control.com



**Physio-Control
Headquarters**
Redmond, WA 98052
www.physio-control.com

Customer Support
P. O. Box 97006
Redmond, WA 98073
Toll free 800 442 1142
Fax 800 426 8049

Physio-Control Canada
Physio-Control Canada Sales, Ltd.
7111 Syntex Drive, 3rd Floor
Mississauga, ON
L5N 8C3
Canada
Toll free 800 895 5896
Fax 866 430 6115



JL Industries, 4450 W 78th St Circle, Bloomington MN 55435
